

Gastric Calibration Tubing



Bariatric
Surgery




Elemental
Healthcare

Advancing Innovation Together

LogiTube™ Gastric Calibration Tubing

For over 30 years Surgical Innovations have been at the forefront of Minimally Invasive Surgery.

LogiTube™ is designed to enable the calibration of the gastric pouch during Sleeve Gastrectomy & Mini Gastric Bypass, or revision surgery for these procedures.



A COST SAVING
OF UP TO
19%¹


1. Data on file



LogiTube™ is equipped with **universal connectivity** to accommodate a variety of syringes, for leak testing, and hospital suction, with finger suction control.



Incremental cm markings allow for precise placement as communicated between the Surgeon and the Anaesthetist.



Distally, six holes are provided, as well as a soft radiopaque tip with a distal suction hole, enabling atraumatic aspiration of stomach contents, or residual fluid, following leak testing.

As healthcare professionals globally try to balance the pent-up demand for surgical interventions, with making healthcare provision economically sustainable, our LogiTube™ range can also deliver a cost saving of up to 19%¹.

LogiTube™ Ordering Information


Product Code	Product Description	Outer Dimensions	Wall Thickness	Total Length	Hole Diameter
34TUBER	34 Fr/Ch Calibration Tubing Red	11.3mm	1.7mm	1,200mm	4.8mm
36TUBEB	36 Fr/Ch Calibration Tubing Blue	12.0mm	2.0mm	1,200mm	4.8mm
38TUBEG	38 Fr/Ch Calibration Tubing Green	12.4mm	2.0mm	1,200mm	4.8mm
40TUBEY	40 Fr/Ch Calibration Tubing Yellow	13.2mm	2.5mm	1,200mm	4.8mm

Leading
the
Evolution
in sustainable surgery

To find out more
information or to request
a demonstration,
visit our website
[www.elementalhealthcare.co.uk/
bariatric-surgery/](http://www.elementalhealthcare.co.uk/bariatric-surgery/)



 Clayton Wood House, 6 Clayton Wood Bank,
Leeds, LS16 6QZ, UK

 +44 (0)844 412 0020

 info@elementalhealthcare.co.uk

 @SurgicalInnovations

 @ElementalHCare

 @ElementalHealthcare

LT-02-23 Version 1